Test Case Planning

Group members	Govan Henry
	Katherine Leon
	Hudson Finocchio
	Jason Gonzalez
	Ty Battle
	Stephen Russell

1. Use Case 2 - Alerts for Medication


Test Case ID	TC01 -
Test Objective	The user correctly sets up a reminder alert
Preconditions	The user must be signed in, have internet connectivity stable, and click on "Reminder"
Test steps	1. Click Add Reminder Alert
	 Fill in the information required for Reminder Alert The user will enter all required fields correctly Hit submit form
Input Values	Medication Name: perc30 (drop down) Reminder time: 3:00 (drop down) Dial Choice AM/PM: PM Start Date: 11/15/2024 (calendar sub-menu) Notification Frequency: Once an hour (drop down) Click the Submit Form button
Expected	Medication will appear below the Add/Update Reminder Notification in the Existing
results	Notifications box Will detail all input information entered and the user's phone number to receive an SMS alert with a delete notification button below it
	* * * * * * * * * * * * * * * * * * *

Test Case ID	TC02 -
Test Objective	The user leaves the AM/PM dial empty for the reminder
Preconditions	The user must be signed in, have internet connectivity stable, and click on "Reminder"
Test steps	User will enter information into each section
	2. The user will not choose the radial dial option for AM/PM
	3. All other sections are filled
	4. The user will click on the submit form
	5. A message should populate to ask to "Please select one of these options." near dial options
Input Values	All values except one field blank Medication Name: Perc30 Reminder time: 1: 00 (drop down) Dial Choice AM/PM: LEFT BLANK Start Date: 11/15/2024 (calendar sub-menu) Notification Frequency: Every 15 mins Click Submit Form
Expected	An error message will display "Please select one of these options." near the AM PM dials
results	● ●
	Add/Update Reminder Notification Medication Name percol Reminder Times 190 - AMN Please select one of these options Start Plater 1/15/2884 Notification Preparety: Every 15 mes Every 15 mes Percol Times 3.00 PM Start Dater 2004-11-15

Test Case ID	TC03 -	
Test Objective	The user leaves the Medication Name blank	
Preconditions	The user must be signed in, have internet connectivity stable, and click on "Reminder"	
Test steps	1. User enters information into the add/ update reminder notification box	
	2. The user will not choose an option for medication name from drop-down	
	3. All other fields will be entered with the desired information	

Test Case ID	TC03 -
	4. Error message displays "Please select an item in the list"
Input Values	Text input: Medication Name: LEFT BLANK Reminder time: 1: 00(drop down) Dial Choice AM/PM: AM Start Date: 11/15/2024(calendar sub-menu) Notification Frequency: Every 15 mins Click Submit Form
Expected results	Displays error message "Please select an item in the list." below Medication Name input

Test Case ID	TC04 -
Test Objective	The user deletes reminders correctly
Preconditions	The user must be signed in, have internet connectivity stable, and click on "Reminder"
Test steps	1. The user scrolls to find medication to delete from the existing medication box
	2. Click on the Delete Reminder button on the desired medication (perc30)
	3. The message populates "Are you sure you want to delete this medication?"
	Clicks on the Yes button
	4. Reminder is deleted
	Existing Notifications
Input Values	Medication: perc30
	Time: 3:00 PM
	Start Date: 2024-11-15
	Frequency: hourly
	Phone: +15403951893
	DELETE NOTIFICATION(Clicked)
	Yes option is clicked on when a message populates



Test Case ID	TC05 -
Test Objective	The user leaves the Start Date blank
Preconditions	The user must be signed in, have internet connectivity stable, and click on "Reminder"
Test steps	1. User will enter information into each section
	2. The user will not choose a Start Date through the calendar button/type in the date
	3. All other sections are filled
	4. The user will click on the submit form
	6. A message should populate to ask to "Please fill out this field" near the Start Date
Input Values	Medication Name: Perc30 Reminder time: 1:00 (drop down) Dial Choice AM/PM: AM Start Date: LEFT BLANK(calendar sub-menu) Notification Frequency: Every 15 mins Click Submit Form
Expected result	Displays message "Please fill out this field." near the Start Date option

